

Infinium® Omni2.5-8 v1.3 BeadChip

Next-generation GWAS content for genotyping and CNV analysis.

Overview

The Infinium Omni2.5-8 v1.3 BeadChip delivers comprehensive coverage of both common and rare SNP content from the 1000 Genomes Project (1kGP; MAF > 2.5%),¹ designed to be maximally informative for diverse world populations. Using the proven HiScan® or iScan® System, these 8-sample BeadChips offer high throughput and optimized tag SNP content, including full support of copy number variation (CNV) applications. The Infinium Omni2.5-8 v1.3 BeadChip is a powerful entry point into the Omni Roadmap, which provides researchers with step-wise, flexible access to 5 million variants per sample. Convenient kit packaging, a streamlined PCR-free protocol, and integrated analysis software are included to provide a comprehensive DNA analysis solution. This product is customizable with up to 200,000 attempted beadtypes with the Infinium Omni2.5-8+ v1.3 version of the BeadChip.



Figure 1: Infinium Omni2.5-8 v1.3 BeadChip—The 8-sample Infinium Omni2.5-8 v1.3 BeadChip supports rapid cost-effective studies with coverage of the latest common and rare variants from the 1000 Genomes Project, down to 2.5% minor allele frequency (MAF).

Table 1: Infinium Omni2.5-8 v1.3 BeadChip Product Information

Feature	Description		
Total Number of Markers	2,372,784		
Capacity for Custom Markers	200,000		
Number of Samples per BeadChip	8		
DNA Input Requirement	200 ng		
Assay	Infinium LCG		
Instrument Support	iScan or HiScan		
Sample Throughput ^a	~1067 samples / week		
Scan Time per Sample	iScan ^a 11.3 min	HiScan 6.5 min	
LD Coverage (<i>r</i> ² ≥ 0.8)			
Population ^b	1kGP ^c MAF > 1%	1kGP ^c MAF > 2.5%	1kGP ^c MAF > 5%
AFR	0.54	0.63	0.69
AMR	0.73	0.82	0.86
EAS	0.79	0.86	0.89
EUR	0.79	0.86	0.89
SAS	0.74	0.83	0.87
Data Performance	Value ^d	Product Specification	
Call Rate	99.8%	> 99% avg.	
Reproducibility	99.99%	> 99.9%	
Log R Deviation	0.11	< 0.30 ^d	
Spacing	Mean	Median	90 th % ^e
Spacing (kb)	1.23	0.65	2.83
<p>a. Estimate assumes 2 iScan systems, one AutoLoader 2.x, 2 Tecan robots, and a 5-day work week.</p> <p>b. See www.1000genomes.org/category/frequently-asked-questions/population</p> <p>c. Compared against the Phase 3 1000 Genomes Project (1kGP) data release (2014).¹</p> <p>d. Values are derived from genotyping 455 HapMap reference samples.</p> <p>e. Value expected for typical projects using standard Illumina protocols. Tumor samples and samples prepared by methods other than standard Illumina protocols are excluded.</p> <p>Abbreviations: 1kGP, 1000 Genomes Project; MAF, minor allele frequency; LD, Linkage Disequilibrium.</p>			

Table 2: Infinium Omni2.5-8 v1.3 BeadChip Marker Information

Marker Categories	Number of Markers ^a		
In RefSeq ^b Genes	1,064,781 (1,272,710 ^c)		
In RefSeq Exons	141,984		
In RefSeq Promoter Regions	50,214		
In ADME Genes	20,848 (26,663 ^c)		
In ADME Exons	2995		
MHC (Extended MHC ^d)	11,116 (18,083)		
Overlap with Genes in COSMIC ^e	966,716		
Overlap with Genes in Gene Ontology ^f	251,141		
Nonsense Markers	375		
Missense Markers	31,986		
Synonymous Markers	32,421		
Silent Markers	77,193		
Mitochondrial Markers	180		
Indels	39		
Sex Chromosomes	X 53,258	Y 2091	Par Loci 2558

a. Compared against the June 2011 1kGP data release.

b. RefSeq - NCBI Reference Sequence Database.²

c. Within 10 kb.

d. MHC is a ~ 4 Mb region; extended MHC is a ~ 8 Mb region.

e. Catalog of somatic mutations in cancer.³

f. Gene Ontology Consortium.⁴

Abbreviations: ADME, absorption, distribution, metabolism, and excretion; MHC, Major Histocompatibility Complex.

Ordering Information

Infinium Omni2.5-8 v1.3 Kit	Catalog No.
16 samples	20001112
48 samples	20001113
96 samples	20001114
384 samples	20001115
Infinium Omni2.5-8+ v1.3 Kit*	Catalog No.
16 samples	20001116
48 samples	20001117
96 samples	20001118
384 samples	20001119
*Enabled for additional custom content	

Learn More

To learn more about the Infinium Omni2.5-8 v1.3 BeadChip and other Illumina genotyping products and services, visit:

www.illumina.com/applications/genotyping.html.

References

1. www.1000genomes.org Accessed 18 April 2014.
2. www.ncbi.nlm.nih.gov/refseq Accessed 6 July 2015.
3. <http://cancer.sanger.ac.uk/cosmic> Accessed 6 July 2015.
4. www.geneontology.org Accessed 6 July 2015.

Illumina • 1.800.809.4566 toll-free (US) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2016 Illumina, Inc. All rights reserved. Illumina, Infinium, iScan, HiScan, and the pumpkin orange color are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. Pub. No. 370-2011-008 Current as of 30 March 2016

